

# Curtin Health Innovation Research Institute – Biosciences

## Flow Cytometry Analyser Guide

### Attune (Acoustic focusing)



#### Attune Flow Cytometer Specifications:

- 2 laser system: Blue / Red Configuration
- Measures 6 fluorochromes simultaneously
- Low fluid consumption (approx. 1L/day)
- Free software
- Sample up to 20,000 events/second
- Directly give cell counts/ml (without beads)
- User friendly software & auto-compensation

## Curtin Health Innovation Research Institute - Biosciences

### Attune

#### Filter Configuration Guide

Channel Name	Standard Filters	Commonly Used Fluorophores
<b>488 nm Blue Laser</b>		
SSC	488/10	Side scatter
[BL1]	530/30	FITC, BB515, AF488, CFSE, GFP, YFP, Live/Dead Fixable Green
[BL2]	574/26	PE
[BL3]	690/50	PerCP, PerCP-Cy5.5, PE-AF700, Live/Dead Fixable
[BL4]	780/60	PE-Cy7, PE-AF750,
<b>638 nm Red Laser</b>		
[RL1]	660/20	APC
[RL2]	780/60	APC-Cy7, APC-H7

**Note:** This is not an exhaustive list of fluorophores, please contact CHIRI Flow Cytometry Staff if you need assistance checking the details of your panel.